



 **VISION**[®]
BIOTECHNOLOGY

“IN QUALITY, WE TRUST”

www.visionbiotechnology.com

ABOUT US

Vision Biyoteknoloji was established in 2012 as Cesna Biotechnology in Kadıköy, Istanbul. Its journey started as a distributor of famous brands including 3M, Merck and Thermo and continued with R&D activities under the brand name Vision Biyoteknoloji in 2015. After developing “Vision GMO Screening Kit” as the first commercial product in 2016, the company very fast became a global Turkish brand, offering services with 240 experts and 160 approved products in 43 countries and 5 continents across the globe, under the motto “On the way of becoming a World brand.”

VISION

Under the motto, “Just Dream, Desire, and Realize”, it is our vision to extent our production and R&D expertise in life sciences and biotechnology to the field of diagnostics. It is our goal to become a globally acknowledged, international brand in diagnostics.

MISSION

Raising awareness of the importance of diagnosis for treatment with the help of our experienced team using the latest technology in our field.

Developing projects and products in the diagnostic field and reaching all parts of the world regardless of the income rate.



REAL-TIME qPCR TEST KITS

Real-time qPCR is a well-established technique that can amplify and quantify or detect the target DNA simultaneously, providing fast and reliable results in diagnostics.

Vision Biyoteknoloji offers a wide range of qPCR test kits for both qualitative and quantitative assessment of various diseases.

INFECTIOUS DISEASES

PRODUCT	QUANTITY	REF
WOUND PANEL KIT	25/50/100 rxn	VSNMP109
WOUND PANEL WITH ANTIBIOTICS RESISTANCE PANEL	25/50/100 rxn	VSNMP108
HEPATITIS B DETECTION TEST KIT	25/50/100 rxn	VSNP001
HEPATITIS C DETECTION TEST KIT	25/50/100 rxn	VSNP002
HIV-1 DETECTION TEST KIT	25/50/100 rxn	VSNP003
HIV-II DETECTION TEST KIT	25/50/100 rxn	VSNHIV02
HBV-HCV DETECTION TEST KIT	25/50/100 rxn	VSNP005
HBV-HCV-HIV-1 MULTIPLEX DETECTION TEST KIT	25/50/100 rxn	VSNP006
HSV-1 DETECTION TEST KIT	25/50/100 rxn	VSNP007
HSV-2 DETECTION TEST KIT	25/50/100 rxn	VSNP008
HSV1-HSV2 DETECTION TEST KIT	25/50/100 rxn	VSNP009
CMV DETECTION TEST KIT	25/50/100 rxn	VSNP010
HEPATITIS D DETECTION KIT	25/50/100 rxn	VSNP011
BKV QUANTIFICATION KIT	25/50/100 rxn	VSNMP002
HHV-6 QUANTIFICATION KIT	25/50/100 rxn	VSNMP005
MTBC DETECTION KIT	25/50/100 rxn	VSNMP008
EBV QUANTIFICATION KIT	25/50/100 rxn	VSNMP009
JCV QUANTIFICATION KIT	25/50/100 rxn	VSNMP011
PARVOVIRUS B19 QUANTIFICATION KIT	25/50/100 rxn	VSNMP012
VZV DETECTION KIT	25/50/100 rxn	VSNMP013
BACTERIAL MENINGITIS SCREENING KIT	25/50/100 rxn	VSNMP014
HELICOBACTER PYLORI DETECTION KIT	25/50/100 rxn	VSNMP019
NEONATAL MENINGITIS SCREENING KIT	25/50/100 rxn	VSNMP020
SEPSIS PATHOGENS SCREENING KIT	25/50/100 rxn	VSNMP021
JCV-BKV DETECTION KIT	25/50/100 rxn	VSNMP031
VESICULAR SCREENING KIT	25/50/100 rxn	VSNMP034
DIPHThERIA DETECTION KIT	25/50/100 rxn	VSNP013
PARASITE INFECTIONS SCREENING KIT	25/50/100 rxn	VSNPI001
MENINGITIS/ENCEPHALIS SCREENING KIT	25/50/100 rxn	VSNMNP01
CANDIDA PCR KIT	25/50/100 rxn	VSNCD001
ENTERIC PCR KIT	25/50/100 rxn	VSNRD016
ENTERIC/ADENO PCR KIT	25/50/100 rxn	VSNRD017
CAMPYLOBACTER PCR KIT	25/50/100 rxn	VSNRD018
CRYPTOCOCCUS NEOFORMANS/CANDIDA PCR KIT	25/50/100 rxn	VSNRD019
CRYPTOCOCCUS SPP DETECTION KIT	25/50/100 rxn	VSNP015

PRODUCT	QUANTITY	REF
CLOSTRIDIUM DIFFICILE DETECTION KIT	25/50/100 rxn	VSNMP038
GASTROENTERITIS SCREENING KIT	25/50/100 rxn	VSNMP040
HCV GENOTYPING KIT	25/50/100 rxn	VSNHCVG01
VIRAL MENINGITIS PANEL KIT	25/50/100 rxn	VSNMP046
VZV (HHV-3) QUANTIFICATION KIT	25/50/100 rxn	VSNMP052
RASH-FEVER SCREENING KIT	25/50/100 rxn	VSNMP053
ENTEROVIRUS DETECTION KIT	25/50/100 rxn	VSNMP054
LEPTOSPIRA SPP. DETECTION KIT	25/50/100 rxn	VSNMP059
CANDIDA ALBICANS DETECTION KIT	25/50/100 rxn	VSNMP061
HHV-8 DETECTION KIT	25/50/100 rxn	VSNMP062
HHV6-HHV7 DETECTION KIT	25/50/100 rxn	VSNMP063
BRUCELLA SPP. DETECTION KIT	25/50/100 rxn	VSNMP066
SMALLPOX DETECTION KIT	25/50/100 rxn	VSNMP067
RUBELLA VIRUS DETECTION KIT	25/50/100 rxn	VSNMP073
THEILERIA SPP. DETECTION KIT	25/50/100 rxn	VSNMP079
BABESIA SPP. DETECTION KIT	25/50/100 rxn	VSNMP080
ANAPLASMA SPP. DETECTION KIT	25/50/100 rxn	VSNMP081
ACINETOBACTER BAUMANNII DETECTION KIT	25/50/100 rxn	VSNMP083
AKKERMANSIA MUCINIPHILA DETECTION KIT	25/50/100 rxn	VSNMP084
BACTEROIDES DETECTION KIT	25/50/100 rxn	VSNMP085
CALDISERICA DETECTION KIT	25/50/100 rxn	VSNMP086
DIALISTER DETECTION KIT	25/50/100 rxn	VSNMP088
URINARY TRACT INFECTIONS KIT	25/50/100 rxn	VSNUT001
WOMEN HEALTH SCREENING KIT	25/50/100 rxn	VSNWHP21
ANTIBIOTICS RESISTANCE SCREENING KIT V1	25/50/100 rxn	VSNAR001
ANTIBIOTICS RESISTANCE SCREENING KIT V2	25/50/100 rxn	VSNAR002
ANTIBIOTICS RESISTANCE SCREENING KIT V3	25/50/100 rxn	VSNAR004
VIBRIO CHOLERAE DETECTION KIT	25/50/100 rxn	VSNVC001
ANTIBIOTIC RESISTANCE/MYCOPLASMA GENITALIUM SCREENING KIT	25/50/100 rxn	VSNAR003
HEPATITIS A-B-C DETECTION KIT	25/50/100 rxn	VSNP014
ACANTHAMOEBA DETECTION KIT	25/50/100 rxn	VSNP030
VANA/VANB ENTEROCOCCUS PCR KIT	25/50/100 rxn	VSNVANAB
BARTONELLA HENSELAE DETECTION KIT	25/50/100 rxn	VSNP022
BRUCELLA TYPING DETECTION KIT	25/50/100 rxn	VSNP023
EBV DETECTION KIT	25/50/100 rxn	VSNMP009-1
EHEC DETECTION KIT	25/50/100 rxn	VSNP026
HELICOBACTER PYLORI CLARITHROMYCIN DETECTION KIT	25/50/100 rxn	VSNMP019-1

PRODUCT	QUANTITY	REF
BACILLUS ANTHRACIS DETECTION KIT	25/50/100 rxn	VSNBA001
YERSINIA PESTIS DETECTION KIT	25/50/100 rxn	VSNYP001
ENTEROCOCCUS SPP. DETECTION KIT	25/50/100 rxn	VSNENT001
WOUND&NAIL SCREENING KIT	25/50/100 rxn	VSNWNP01
BONE&JOINT INFECTIONS SCREENING KIT	25/50/100 rxn	VSNBJ001
HEPATITIS A DETECTION KIT	25/50/100 rxn	VSNP033
HEPATITIS G DETECTION TEST KIT	25/50/100 rxn	VSNP028
URINARY TRACT INFECTIONS KIT V2	25/50/100 rxn	VSNUT002

GENETIC DISEASES

PRODUCT	QUANTITY	REF
FMF MULTIPLEX REAL-TIME PCR TEST KIT	25/50/100 rxn	VISNG001
FMF MUTATION ANALYSIS KIT	25/50/100 rxn	VSNMP004
K-28B POLYMORPHISM DETECTION KIT	25/50/100 rxn	VISNG003
PAI-1 4G/5G DETECTION KIT	25/50/100 rxn	VSNMP003
HLAB27 DETECTION KIT	25/50/100 rxn	VSNMP006
FXIII DETECTION KIT	25/50/100 rxn	VSNMP010
PROTHROMBIN DETECTION KIT	25/50/100 rxn	VSNMP022
MTHFR C677T DETECTION KIT	25/50/100 rxn	VSNMP115
FVL DETECTION KIT	25/50/100 rxn	VSNMP026
MTHFR A1298C DETECTION KIT	25/50/100 rxn	VSNMP028
HLA-B*57:01 DETECTION KIT	25/50/100 rxn	VSNMP030
THROMBOPHILIA DETECTION KIT	25/50/100 rxn	VSNMP037
THROMBOPHILIA MULTIPLEX DETECTION KIT V2	25/50/100 rxn	VSNMP037-2
FVC DETECTION KIT	25/50/100 rxn	VSNMP047
FII-FVL DETECTION KIT	25/50/100 rxn	VSNMP050
FV H1299R DETECTION KIT	25/50/100 rxn	VSNMP051
SMA REAL TIME PCR KIT	25/50/100 rxn	VSNSMA01

RESPIRATORY PANEL

PRODUCT	QUANTITY	REF
COVID19 EASYPREP FAST REAL-TIME PCR TEST KIT	25/50/100 rxn	V2004F190
SARS-COV-2 MUTATION DETECTION KIT	25/50/100 rxn	V2102HG24
SARS-COV-2/FLU DETECTION KIT	25/50/100 rxn	VSNRD001
RESPIRATORY PATHOGENS PANEL(RPP) MINI	25/50/100 rxn	VSNMP117
RESPIRATORY PATHOGENS PANEL(RPP) MIDI	25/50/100 rxn	VSNMP118
RESPIRATORY PATHOGENS PANEL(RPP) MAXI	25/50/100 rxn	VSNMP119
HUMAN METAPNEUMOVIRUS DETECTION KIT	25/50/100 rxn	VSNR001
WHITE LUNG SYNDROME MINI PANEL PCR TEST KIT	25/50/100 rxn	VSNRD015
WHITE LUNG SYNDROME MIDDLE PANEL PCR TEST KIT	25/50/100 rxn	VSNRD022
MYCOPLASMA PNEUMONIAE DETECTION KIT	25/50/100 rxn	VSNR002
LEGIONELLA PNEUMOPHILA DETECTION KIT	25/50/100 rxn	VSNR003
KLEBSIELLA PNEUMONIAE DETECTION KIT	25/50/100 rxn	VSNR004
BORDETELLA PERTUSSIS DETECTION KIT	25/50/100 rxn	VSNR005
CHLAMYDIA PNEUMONIAE DETECTION KIT	25/50/100 rxn	VSNR006
PNEUMOCYSTIS JIROVECI DETECTION KIT	25/50/100 rxn	VSNR007
RESPIRATORY BACTERIAL BASIC SCREENING KIT	25/50/100 rxn	VSNMP048-B
INFLUENZA A- H5N1 DETECTION KIT	25/50/100 rxn	VSNINF001
LEGIONELLA MULTIPLEX DETECTION KIT	25/50/100 rxn	VSNLP001
MORAXELLA CATARRHALIS DETECTION KIT	25/50/100 rxn	VSNR008
RESPIRATORY SYNCYTIAL VIRUS (RSV) – A DETECTION KIT	25/50/100 rxn	VSNR009
RESPIRATORY SYNCYTIAL VIRUS (RSV) – B DETECTION KIT	25/50/100 rxn	VSNR010
RESPIRATORY SYNCYTIAL VIRUS (RSV) – A,B DETECTION KIT	25/50/100 rxn	VSNR011
INFLUENZA A DETECTION KIT	25/50/100 rxn	VSNR012
INFLUENZA B DETECTION KIT	25/50/100 rxn	VSNR013
INFLUENZA-A,B DETECTION KIT	25/50/100 rxn	VSNR014
ADENOVIRUS DETECTION KIT	25/50/100 rxn	VSNR015
BOCAVIRUS DETECTION KIT	25/50/100 rxn	VSNR016
RHINOVIRUS (A, B, C) DETECTION KIT	25/50/100 rxn	VSNR017
PARAINFLUENZA (1, 2, 3, 4) DETECTION KIT	25/50/100 rxn	VSNR018
HUMAN PARECHOVIRUS DETECTION KIT	25/50/100 rxn	VSNR019
ASPERGILLUS DETECTION KIT	25/50/100 rxn	VSNMP007
H1N1 DETECTION KIT	25/50/100 rxn	VSNMP025
INFLUENZA A-B-C SCREENING KIT	25/50/100 rxn	VSNMP033
INFLUENZA TYPING KIT	25/50/100 rxn	VSNMP041

PRODUCT**QUANTITY****REF**

RSV-FLU DETECTION KIT	25/50/100 rxn	VSNMP042
ADEP SCREENING KIT	25/50/100 rxn	VSNMP044
ADENOVIRUS QUANTIFICATION KIT	25/50/100 rxn	VSNMP045
RESPIRATORY VIRAL BASIC SCREENING KIT	25/50/100 rxn	VSNMP048
B. PERTUSSIS/PARAPERTUSSIS GENOTYPING KIT	25/50/100 rxn	VSNMP055
STREPTOCOCCUS PNEUMONIAE DETECTION KIT	25/50/100 rxn	VSNMP058
STAPHYLOCOCCUS AUREUS DETECTION KIT	25/50/100 rxn	VSNMP060
MERS-COV DETECTION KIT	25/50/100 rxn	VSNMP071
STREPTOCOCCUS PYOGENES DETECTION KIT	25/50/100 rxn	VSNMP078
HAEMOPHILUS INFLUENZA TYPING (A/B/C/D/E/F) DETECTION KIT	25/50/100 rxn	VSNMP082
STREPTOCOCCUS A/C/G PCR TEST KIT	25/50/100 rxn	VSNSTR01
SWINE FLU DETECTION KIT	25/50/100 rxn	VSNSF001
STREP C/G DETECTION KIT	25/50/100 rxn	VSNMP078-1

TROPICAL AND OTHER VIRUSES, BACTERIA AND PARASITES**PRODUCT****QUANTITY****REF**

MONKEYPOX DETECTION KIT	25/50/100 rxn	VISMPX001
LYME DISEASE DETECTION KIT	25/50/100 rxn	VSNTD001
DENGUE VIRUS DETECTION KIT	25/50/100 rxn	VSNTD002
MALARIA DETECTION KIT	25/50/100 rxn	VSNTD003
WEST NILE DETECTION KIT	25/50/100 rxn	VSNTD004
ZIKA VIRUS DETECTION KIT	25/50/100 rxn	VSNTD005
CHIKUNGUNYA VIRUS DETECTION KIT	25/50/100 rxn	VSNTD006
WEST NILE VIRUS QUANTIFICATION KIT	25/50/100 rxn	VSNMP001
DENGUE VIRUS GENOTYPING KIT	25/50/100 rxn	VSNTD018
TROPICAL FEVER SCREENING KIT	25/50/100 rxn	VSNMP032
MALARIA GENOTYPING KIT	25/50/100 rxn	VSNMP035
CCHFV QUANTIFICATION KIT	25/50/100 rxn	VSNMP039
YELLOW FEVER DETECTION KIT	25/50/100 rxn	VSNMP069
ZIKA VIRUS DETECTION KIT	25/50/100 rxn	VSNMP070
EBOLA VIRUS DETECTION KIT	25/50/100 rxn	VSNMP072
CCHFV DETECTION KIT	25/50/100 rxn	VSNMP039-1
HELMINTHS REAL-TIME PCR KIT	25/50/100 rxn	VSNH001

STDs

PRODUCT	QUANTITY	REF
GARDNERELLA VAGINALIS DETECTION KIT	25/50/100 rxn	VSNS001
TRICHOMONAS VAGINALIS DETECTION KIT	25/50/100 rxn	VSNS002
TREPONEMA PALLIDUM DETECTION KIT	25/50/100 rxn	VSNS003
MYCOPLASMA HOMINIS DETECTION KIT	25/50/100 rxn	VSNS004
MYCOPLASMA GENITALIUM DETECTION KIT	25/50/100 rxn	VSNS005
NEISSERIA GONORRHOEAE DETECTION KIT	25/50/100 rxn	VSNS006
HPV SCREENING KIT V1	25/50/100 rxn	VSNP004-E
HPV (16,18) HIGH RISK DETECTION KIT	25/50/100 rxn	VSNP004-B
HPV (16,18,31) HIGH RISK DETECTION KIT	25/50/100 rxn	VSNP004-C
HPV GENOTYPING REAL-TIME PCR KIT	25/50/100 rxn	VSNP004-D
HPV SCREENING KIT V2	25/50/100 rxn	VSNP004-F
HPV FULL GENOTYPING KIT	25/50/100 rxn	VSNP004-G
CHLAMYDIA TRACHOMATIS DETECTION KIT	25/50/100 rxn	VSNR020
HPV 6-11-16-18 LR-HR DETECTION KIT	25/50/100 rxn	VSNMP065
STDS SCREENING KIT	25/50/100 rxn	VSNMP015
UREAPLASMA UREALYTICUM DETECTION KIT	25/50/100 rxn	VSNMP049
BVAB (GV/CA/TV) PCR KIT	25/50/100 rxn	VSNMP113
UREAPLASMA PARVUM DETECTION KIT	25/50/100 rxn	VSNP017
UREAPLASMA UREALYTICUM/PARVUM DETECTION KIT	25/50/100 rxn	VSNP018
SEXUALLY TRANSMITTED INFECTIONS KIT 28T	25/50/100 rxn	VSNSTD003
SEXUALLY TRANSMITTED INFECTIONS KIT 14T	25/50/100 rxn	VSNSTD004
STI SCREENING KIT	25/50/100 rxn	VSNSTIP10

LIFE SCIENCE

PRODUCT	QUANTITY	REF
VISOL ISOLATION REAGENT	100 rxn	VSNMP074
VIS-MIX ONE STEP RT-QPCR PROBE MASTER MIX 2X	100 rxn	VSNMP075
ONE TAQ HOT START 2X MASTER MIX	100 rxn	VSNMP076
CDNA SYNTHESIS KIT	50 rxn	VSNLS001

MICROBIOTA

PRODUCT	QUANTITY	REF
LACHNOSPIRA DETECTION KIT	25/50/100 rxn	VSNMP095
LACHNOSPIRACEAE DETECTION KIT	25/50/100 rxn	VSNMP096
LACTOBACILLUS ACIDOPHILUS DETECTION KIT	25/50/100 rxn	VSNMP097
LACTOBACILLUS BREVIS DETECTION KIT	25/50/100 rxn	VSNMP098
LACTOBACILLUS FERMENTUM DETECTION KIT	25/50/100 rxn	VSNMP099
LACTOBACILLUS PLANTARUM DETECTION KIT	25/50/100 rxn	VSNMP100
LACTOBACILLUS RHAMNOSUS DETECTION KIT	25/50/100 rxn	VSNMP101
PREVOTELLACEAE DETECTION KIT	25/50/100 rxn	VSNMP102
PROTEOBACTERIA DETECTION KIT	25/50/100 rxn	VSNMP103
SALMONELLA ENTERITIDIS DETECTION KIT	25/50/100 rxn	VSNMP104
SALMONELLA TYPHIMURIUM DETECTION KIT	25/50/100 rxn	VSNMP105
SUBDOLIGRANULUM DETECTION KIT	25/50/100 rxn	VSNMP106
SUCCINIVIBRIO DETECTION KIT	25/50/100 rxn	VSNMP107
TENERICUTES DETECTION KIT	25/50/100 rxn	VSNMP108

BYPRODUCTS

PRODUCT	QUANTITY	REF
VIS-EXPRESS	100	VL301711
VIRAL TRANSPORT MEDIUM	100	VVTM001
SMEAR SAMPLING SET	20	VSME

FOOD PANEL

PRODUCT	QUANTITY	REF
GMO SCREENING KIT	100 rxn	VSNF001
PLANT SPECIFIC KIT	100 rxn	VSNF002
SOY-COTTON SCREENING KIT	100 rxn	VSNF003
EURL FORMAT SOY EVENT MON87701 OLIGONUCLEOTIDE SET	100 rxn	VSNF004

FOOD PANEL

PRODUCT

QUANTITY

REF

EURL FORMAT SOY EVENT MON87708 OLIGONUCLEOTIDE SET	100 rxn	VSNF005
EURL FORMAT SOY EVENT MON87751 OLIGONUCLEOTIDE SET	100 rxn	VSNF006
EURL FORMAT SOY EVENT MON87769 OLIGONUCLEOTIDE SET	100 rxn	VSNF007
EURL FORMAT SOY EVENT CV127 OLIGONUCLEOTIDE SET	100 rxn	VSNF008
EURL FORMAT SOY EVENT DP-305423 OLIGONUCLEOTIDE SET	100 rxn	VSNF009
EURL FORMAT SOY EVENT DAS-44406-6 OLIGONUCLEOTIDE SET	100 rxn	VSNF010
EURL FORMAT SOY EVENT DAS-68416-4 OLIGONUCLEOTIDE SET	100 rxn	VSNF011
EURL FORMAT SOY EVENT DAS-81419-2 OLIGONUCLEOTIDE SET	100 rxn	VSNF012
EURL FORMAT COTTON EVENT 281-24-236 OLIGONUCLEOTIDE SET	100 rxn	VSNF013
EURL FORMAT COTTON EVENT 3006-210-23 OLIGONUCLEOTIDE SET	100 rxn	VSNF014
EURL FORMAT COTTON EVENT GHB614 OLIGONUCLEOTIDE SET	100 rxn	VSNF015
EURL FORMAT COTTON EVENT GHB811 OLIGONUCLEOTIDE SET	100 rxn	VSNF023
EURL FORMAT MAIZE HMG DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF016
EURL FORMAT MAIZE ADH1 DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF035
EURL FORMAT MAIZE EVENT DAS40278-9 OLIGONUCLEOTIDE SET	100 rxn	VSNF017
SOY EVENT 1 (MON87701, MON87708, MON87769) SCREENING KIT	100 rxn	VSNF018
SOY EVENT 2 (CV127- DP305423) SCREENING KIT	100 rxn	VSNF019
SOY EVENT 3 (DAS-44406-6 / DAS-68416-4 / DAS-81419-2) SCREENING KIT	100 rxn	VSNF020
COTTON EVENT (281-24-236, 3006-210-23, GHB614) SCREENING KIT	100 rxn	VSNF021
EURL FORMAT SOY LECTIN DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF024
EURL FORMAT COTTON ACP1 DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF025
EURL FORMAT COTTON SAH7 DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF036
EURL FORMAT CANOLA EVENT DP-073496-4 OLIGONUCLEOTIDE SET	100 rxn	VSNF034
EURL FORMAT CANOLA CRUA DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF026
EURL FORMAT P35S DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF027
EURL FORMAT PFMV DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF028
EURL FORMAT TNOS DETECTION OLIGONUCLEOTIDE SET	100 rxn	VSNF029
LISTERIA MONOCYTOGENES DETECTION KIT	100 rxn	VSNF031
GLUTEN DETECTION KIT	100 rxn	VSNF033
MAIZE-CANOLA SCREENING KIT	100 rxn	VSNF022
SALMONELLA SPP. DETECTION KIT	100 rxn	VSNMP068
ELUSIMICROBIA DETECTION KIT	100 rxn	VSNMP089
ESCHERICHIA COLI O157:H7 DETECTION KIT	100 rxn	VSNMP090
FAECALIBACTERIUM PRAUSNITZII DETECTION KIT	100 rxn	VSNMP092
FAECALIBACTERIUM DETECTION KIT	100 rxn	VSNMP093



LFA RAPID TEST KITS

Lateral Flow Immunoassay, or rapid test, is a fast cheap, and simple device intended to detect the presence of a target substance.

LFAs are widely used in medical diagnostics in the home, at the point of care, and in the laboratory.

Vision Biyoteknoloji offers a wide range of rapid test kits for the qualitative assessment of various diseases.

COVID-19

PRODUCT	QUANTITY	REF
COVID-19 Ag RAPID TEST KIT	1/25/50/100	V2004Y190
COVID-19 POC Ab RAPID TEST KIT	1/25/50/100	V2004R190
SARS-CoV-2 NEUTRALIZING ANTIBODY RAPID TEST KIT	1/25/50/100	VISV2207AGT061

CARDIAC MARKERS

PRODUCT	QUANTITY	REF
CK-MB RAPID TEST KIT	1/25/50/100	VISD011
cTnI-TROPONIN Ag RAPID TEST KIT	1/25/50/100	VISD004
Myo- cTnI RAPID TEST KIT	1/25/50/100	VISD021
Myo/CK-MB/cTnI RAPID TEST KIT	1/25/50/100	VISD049
MYOGLOBIN Ag RAPID TEST KIT	1/25/50/100	VISD006

RESPIRATORY DISEASE

PRODUCT	QUANTITY	REF
COVID-INF A/B- RSV- ADENO MULTI RAPID TEST KIT	1/25/50/100	VISD060
COVID-INF A/B-RSV MULTI RAPID TEST KIT	1/25/50/100	VISD059
INFLUENZA A/B RAPID TEST KIT	1/25/50/100	VISD007
INFLUENZA/COVID RAPID TEST KIT	1/25/50/100	VISD056
MYCOPLASMA RAPID TEST KIT	1/25/50/100	VISD047
RSV RAPID TEST KIT	1/25/50/100	VISD023
STREP A RAPID TEST KIT	1/25/50/100	VISD024
TUBERCULOSIS RAPID TEST KIT	1/25/50/100	VISD016
WHITE LUNG SYNDROME MINI PANEL RAPID TEST KIT	1/25/50/100	VISD078
WHITE LUNG SYNDROME MIDDLE PANEL RAPID TEST KIT	1/25/50/100	VISD079
INFLUENZA A-B/RSV/ADENO RAPID TEST KIT	1/25/50/100	VISD087
FLU A/B-RSV RAPID TEST KIT	1/25/50/100	VISD088
HMPV (HUMAN METAPNEUMOVIRUS) ANTIGEN RAPID TEST KIT	1/25/50/100	VISD121
COVID-INF A/B- RSV- ADENO- HMPV MULTI RAPID TEST KIT	1/25/50/100	VISD122

INFECTIOUS DISEASES

PRODUCT	QUANTITY	REF
CHIK Ab RAPID TEST KIT	1/25/50/100	VISD031
CMV Ab RAPID TEST KIT	1/25/50/100	VISD028
DENGUE IgG/IgM RAPID TEST KIT	1/25/50/100	VISD025
DENGUE NS1&IgG/IgM RAPID TEST KIT	1/25/50/100	VISD001
DENGUE AG RAPID TEST KIT	1/25/50/100	VISD035
HSV I Ab RAPID TEST KIT	1/25/50/100	VISD038
HSV II Ab RAPID TEST KIT	1/25/50/100	VISD039
HSV I-II Ab RAPID TEST KIT	1/25/50/100	VISD043
LEGIONELLA RAPID TEST KIT	1/25/50/100	VISD036
MALARIA Pf RAPID TEST KIT	1/25/50/100	VISD002
MALARIA Pf/Pan RAPID TEST KIT	1/25/50/100	VISD019
MALARIA Pf/Pv RAPID TEST KIT	1/25/50/100	VISD005
MALARIA Pf/Pv/Pan RAPID TEST KIT	1/25/50/100	VISD045
MULTI INFECTION RAPID TEST KIT	1/25/50/100	VISD062
RUBELLA Ab RAPID TEST KIT	1/25/50/100	VISD027
TORCH IgG/IgM COMBO RAPID TEST KIT	1/25/50/100	VISD048
TOXOPLASMOSIS Ab RAPID TEST KIT	1/25/50/100	VISD026
TYHPOID RAPID TEST KIT	1/25/50/100	VISD037
YELLOW FEVER RAPID TEST KIT	1/25/50/100	VISD063
ZIKA Ab RAPID TEST KIT	1/25/50/100	VISD030

STDs

PRODUCT	QUANTITY	REF
HBsAg RAPID TEST KIT	1/25/50/100	VISD022
HBsAb RAPID TEST KIT	1/25/50/100	VISD008
HCV Ab RAPID TEST KIT	1/25/50/100	VISD017
ANTI-TP (SYPHILIS) RAPID TEST KIT	1/25/50/100	VISD003
CHLAMYDIA Ab RAPID TEST KIT	1/25/50/100	VISD061
CHLAMYDIA Ag RAPID TEST KIT	1/25/50/100	VISD029
GONORRHEA Ag RAPID TEST KIT	1/25/50/100	VISD018
HIV I-II Ab RAPID TEST KIT	1/25/50/100	VISD010
HPV (16-18) HIGH RISK RAPID TEST KIT	1/25/50/100	VISD058
HPV (6-11) LOW RISK RAPID TEST KIT	1/25/50/100	VISD062
HAsAb RAPID TEST KIT	1/25/50/100	VISD082
HEsAb RAPID TEST KIT	1/25/50/100	VISD083

DRUG TESTS

PRODUCT	QUANTITY	REF
ALCOHOL RAPID TEST KIT	1/25/50/100	VISSETG01
AMPHETAMINE RAPID TEST KIT	1/25/50/100	VISAMP01
BARBITURATES RAPID TEST KIT	1/25/50/100	VISBAR01
BENZODIAZEPINE RAPID TEST KIT	1/25/50/100	VISBEN01
BUPRENORPHINE RAPID TEST KIT	1/25/50/100	VISBUP01
COCAINE RAPID TEST KIT	1/25/50/100	VISCOC01
FENTANYL RAPID TEST KIT	1/25/50/100	VISFEN01
K2 SPICE RAPID TEST KIT	1/25/50/100	VISK201
KETAMINE RAPID TEST KIT	1/25/50/100	VISKET01
METHADONE RAPID TEST KIT	1/25/50/100	VISMETHA01
METHAMPHETAMINE RAPID TEST KIT	1/25/50/100	VISMET01
MDMA RAPID TEST KIT	1/25/50/100	VISMDMA01
MORPHINE RAPID TEST KIT	1/25/50/100	VISMOP01
MULTI DRUG RAPID TEST KIT 6	1/25/50/100	VISMD02
MULTI DRUG RAPID TEST KIT 9	1/25/50/100	VISMD03
MULTI DRUG RAPID TEST KIT 12	1/25/50/100	VISMD04
MULTI DRUG RAPID TEST CUP	1/25/50/100	VISMD05
NICOTINE RAPID TEST KIT	1/25/50/100	VISCOT01
OPIUM RAPID TEST KIT	1/25/50/100	VISOPI01
OXYDIHYDROCODEINONE RAPID TEST KIT	1/25/50/100	VISOXY01
PHENCYLIDINE RAPID TEST KIT	1/25/50/100	VISPCP01
TETRAHYDROCANNOBINOL RAPID TEST KIT	1/25/50/100	VISTET01
TRICYCLIC ANTIDEPRESSANTS RAPID TEST KIT	1/25/50/100	VISTCA01
SALIVA ALCOHOL TEST STRIP	1/25/50/100	VISSETG02

TUMOR MARKERS TEST

PRODUCT	QUANTITY	REF
FOB RAPID TEST KIT	1/25/50/100	VISI001

FERTILITY TEST

PRODUCT	QUANTITY	REF
PREGNANCY TEST KIT	1/25/50/100	VISD103
OVULATION RAPID TEST KIT	1/25/50/100	VISD096

DIGESTION TESTS

PRODUCT	QUANTITY	REF
ADENOVIRUS RAPID TEST KIT	1/25/50/100	VISD014
CHOLERA RAPID TEST KIT	1/25/50/100	VISD057
CRYPTOSPORIDIUM RAPID TEST KIT	1/25/50/100	VISD041
CRYPTOSPORIDIUM-GIARDIA-ENTAMOEBA RAPID TEST KIT	1/25/50/100	VISD040
ENTAMOEBA RAPID TEST KIT	1/25/50/100	VISD042
GIARDIA RAPID TEST KIT	1/25/50/100	VISD032
H. PYLORI Ab RAPID TEST KIT	1/25/50/100	VISD013
H. PYLORI Ag RAPID TEST KIT	1/25/50/100	VISD012
NOROVIRUS RAPID TEST KIT	1/25/50/100	VISD033
ROTAVIRUS Ag RAPID TEST KIT	1/25/50/100	VISD015
ROTAVIRUS-ADENOVIRUS Ag RAPID TEST KIT	1/25/50/100	VISD009
TOXIN A/B RAPID TEST KIT	1/25/50/100	VISD034
URINE STRIP RAPID TEST KIT	1/25/50/100	VISD044
C. DIFFICILE TOXIN A/TOXIN B/GDH RAPID TEST KIT	1/25/50/100	VISD075

MILK TESTS

PRODUCT	QUANTITY	REF
BETA-LACTAM ANTIBIOTIC RAPID TEST KIT	96	VISDF001
BETA-LACTAM/TETRACYCLINE ANTIBIOTIC RAPID TEST KIT	96	VISDF002
BETA-LACTAM/TETRACYCLINE/CEFTIOFUR ANTIBIOTIC RAPID TEST KIT	96	VISDF003
BETA-LACTAMS/TETRACYCLINES/CEPHALEXIN ANTIBIOTIC RAPID TEST KIT	96	VISDF004
BETA-LACTAMS/CEPHALEXIN RAPID TEST KIT	96	VISDF005
BETA-LACTAMS/TETRACYCLINES/CEPHALEXIN/MELAMINE RAPID TEST KIT	96	VISDF006
BETA-LACTAMS/TETRACYCLINES/CHLORAMPHENICOL/STREPTOMYCIN RAPID TEST KIT	96	VISDF007
MELAMINE RAPID TEST KIT	96	VISDF008
GENTAMICIN RAPID TEST KIT	96	VISDF009
KANAMYCIN RAPID TEST KIT	96	VISDF010
NEOMYCIN RAPID TEST KIT	96	VISDF011
SPECTINOMYCIN RAPID TEST KIT	96	VISDF012
STREPTOMYCIN RAPID TEST KIT	96	VISDF013
STREPTOMYCIN/DIHYDROSTREPTOMYCIN RAPID TEST KIT	96	VISDF014
STREPTOMYCIN/DIHYDROSTREPTOMYCIN/GENTAMICIN RAPID TEST KIT	96	VISDF015
CHLORAMPHENICOL RAPID TEST KIT	96	VISDF016
STREPTOMYCIN/DIHYDROSTREPTOMYCIN/GENTAMICIN/NEOMYCIN RAPID TEST KIT	96	VISDF017

PRODUCT	QUANTITY	REF
CHLORAMPHENICOL/CHLORAMPHENICOL SUCCINATE RAPID TEST KIT	96	VISDF018
GENTAMICIN/NEOMYCIN/STREPTOMYCIN/SPECTINOMYCIN/KANAMYCIN RAPID TEST KIT	96	VISDF019
SULFONAMIDE RAPID TEST KIT	96	VISDF020
SULFONAMIDE-TRIMETHOPRIM RAPID TEST KIT	96	VISDF021
AFLATOXIN B1 RAPID TEST KIT	96	VISDF022
AFLATOXIN M1 RAPID TEST KIT	96	VISDF023
BETAMETHASONE RAPID TEST KIT	96	VISDF024
DEXAMETHASONE RAPID TEST KIT	96	VISDF025
DEXAMETHASONE/BETAMETHASONE RAPID TEST KIT	96	VISDF026
ERYTHROMYCIN RAPID TEST KIT	96	VISDF027
SPIRAMYCIN RAPID TEST KIT	96	VISDF028
TYLOSIN RAPID TEST KIT	96	VISDF029
ERYTHROMYCIN/SPIRAMYCIN/TYLOSIN RAPID TEST KIT	96	VISDF030
QUINOLONES RAPID TEST KIT	96	VISDF031
SALBUTAMOL RAPID TEST KIT	96	VISDF032
ZEARAITOL RAPID TEST KIT	96	VISDF033
ZEARALENOL-ZEARALENONE RAPID TEST KIT	96	VISDF034
BACITRACIN RAPID TEST KIT	96	VISDF035
CYPROMAZINE RAPID TEST KIT	96	VISDF036
CLAVULANIC ACID RAPID TEST KIT	96	VISDF037
RACTOPAMINE RAPID TEST KIT	96	VISDF038
DEOXYNIVALENOL (DON) RAPID TEST KIT	96	VISDF039
BETA-LACTAMASE RAPID TEST KIT	96	VISDF040
FLORFENICOL RAPID TEST KIT	96	VISDF041
THIAMPHENICOL RAPID TEST KIT	96	VISDF042
FLORFENICOL-THIAMPHENICOL RAPID TEST KIT	96	VISDF043
ALBENDAZOLE RAPID TEST KIT	96	VISDF044
ABAMECTIN RAPID TEST KIT	96	VISDF045
ANALGIN RAPID TEST KIT	96	VISDF046
OFENDAZOLE RAPID TEST KIT	96	VISDF047
BENZOIC ACID RAPID TEST KIT	96	VISDF048
BETA-AGONIST (4IN1) RAPID TEST KIT	96	VISDF049
BETA-AGONIST (6IN1) RAPID TEST KIT	96	VISDF050
BETA-AGONIST (9IN1) RAPID TEST KIT	96	VISDF051
MELOXICAM RAPID TEST KIT	96	VISDF052
DORAMECTIN RAPID TEST KIT	96	VISDF053
FLORFENICOLAMIDE RAPID TEST KIT	96	VISDF054

PRODUCT**QUANTITY****REF**

TILMICOSIN RAPID TEST KIT	96	VISDF055
ACETAMINOPHEN RAPID TEST KIT	96	VISDF056
FIPRONIL RAPID TEST KIT	96	VISDF057
FLUNIXIN MEGLUMINE RAPID TEST KIT	96	VISDF058
LINCOMYCIN RAPID TEST KIT	96	VISDF059
GOAT MILK ADULTERATION TEST KIT	96	VISDF060

VETERINARY TESTS**PRODUCT****QUANTITY****REF**

CANINE DISTEMPER VIRUS (CDV) Ag RAPID TEST KIT	25/50/100	VISDV001
CANINE PARVOVIRUS (CPV) Ag RAPID TEST KIT	25/50/100	VISDV002
CORONAVIRUS (CCV) Ag RAPID TEST KIT	25/50/100	VISDV003
CANINE INFLUENZA VIRUS (CIV) Ag RAPID TEST KIT	25/50/100	VISDV004
CANINE PARVOVIRUS/CORONAVIRUS (CPV/CCV) Ag MULTI RAPID TEST KIT	25/50/100	VISDV005
RABIES Ag RAPID TEST KIT	25/50/100	VISDV006
PREGNANCY RELAXIN (RLN) RAPID TEST KIT	25/50/100	VISDV007
LEISHMANIA CANIS (LSH) Ab RAPID TEST KIT	25/50/100	VISDV008
CANINE BRUCELLA Ab RAPID TEST KIT	25/50/100	VISDV009
CANINE GIARDIA Ag RAPID TEST KIT	25/50/100	VISDV010
CANINE TOXOPLASMA Ab RAPID TEST KIT	25/50/100	VISDV011
FELINE PANLEUCOPENIA (FPV) Ag RAPID TEST KIT	25/50/100	VISDV012
FELINE CORONAVIRUS (FCoV) Ag RAPID TEST KIT	25/50/100	VISDV013
FELINE CALICIVIRUS (FCV) Ag RAPID TEST KIT	25/50/100	VISDV014
FELINE VIRAL RHINOTRACHEITIS (FHV) Ag RAPID TEST KIT	25/50/100	VISDV015
FELINE INFECTIOUS PERITONITIS (FIP) Ag RAPID TEST KIT	25/50/100	VISDV016
FELINE INFECTIOUS PERITONITIS (FIP) Ab RAPID TEST KIT	25/50/100	VISDV017
FELINE IMMUNODEFICIENCY VIRUS (FIV) Ab RAPID TEST KIT	25/50/100	VISDV018
FIV Ab/FeLV Ag MULTI RAPID TEST KIT	25/50/100	VISDV019
FELINE GIARDIA Ag RAPID TEST KIT	25/50/100	VISDV020
FELINE TOXOPLASMA Ab RAPID TEST KIT	25/50/100	VISDV021
FELINE LEUKEMIA VIRUS (FeLV) Ag RAPID TEST KIT	25/50/100	VISDV022

OTHER TESTS

PRODUCT	QUANTITY	REF
TETANUS IGG RAPID TEST KIT	25/50/100	VISD067
VITAMIN D RAPID TEST KIT	25/50/100	VISD065
SALMONELLA S. RAPID TEST KIT	25/50/100	VISD069
FILARIASIS RAPID TEST KIT	25/50/100	VISD080
LEPTOSPIRA RAPID TEST KIT	25/50/100	VISD081



FIA TEST KITS

VISION FIA TESTS has the ability to perform qualitative analysis with high sensitivity and specificity using fluorescent materials.

It saves time thanks to automatic strip recognition and early detection by reading the information stored in the 2D barcode with VISION fluorescent immunoassay analyzer.

PRODUCT**QUANTITY****REF**

cTnl RAPID TEST KIT (FIA)	25/50/100	VFIA-CTNL
hs-cTnl RAPID TEST KIT (FIA)	25/50/100	VFIA-HSCTNL
Myo RAPID TEST KIT	25/50/100	VFIA-MYO
CK-MB RAPID TEST KIT	25/50/100	VFIA-CK
NT-proBNP RAPID TEST KIT	25/50/100	VFIA-NT
D-Dimer RAPID TEST KIT	25/50/100	VFIA-DD
PCT RAPID TEST KIT	25/50/100	VFIA-PCT
IL-6 RAPID TEST KIT (FIA)	25/50/100	VFIA-IL
CRP RAPID TEST KIT (FIA)	25/50/100	VFIA-CRP
SAA RAPID TEST KIT	25/50/100	VFIA-SAA
HbA1c RAPID TEST KIT	25/50/100	VFIA-HBA
25-OH-VD RAPID TEST KIT	25/50/100	VFIA-D
FERRITIN RAPID TEST KIT	25/50/100	VFIA-FRT
VB12 RAPID TEST KIT	25/50/100	VFIA-VB21
FA RAPID TEST KIT (FIA)	25/50/100	VFIA-FA
C-PEPTIDE RAPID TEST KIT	25/50/100	VFIA-CPEP
INSULIN RAPID TEST KIT	25/50/100	VFIA-INS
CEA RAPID TEST KIT	25/50/100	VFIA-CEA
fPSA RAPID TEST KIT	25/50/100	VFIA-FPSA
PSA RAPID TEST KIT	25/50/100	VFIA-PSA
AMH RAPID TEST KIT (FIA)	25/50/100	VFIA-AMH
LH RAPID TEST KIT (FIA)	25/50/100	VFIA-LH
FSH RAPID TEST KIT	25/50/100	VFIA-FSH
PRL RAPID TEST KIT	25/50/100	VFIA-PRL
TES RAPID TEST KIT	25/50/100	VFIA-TES
PROG RAPID TEST KIT	25/50/100	VFIA-PROG
TSH RAPID TEST KIT	25/50/100	VFIA-TSH
T3 RAPID TEST KIT	25/50/100	VFIA-T3
T4 RAPID TEST KIT	25/50/100	VFIA-T4
FT3 RAPID TEST KIT	25/50/100	VFIA-FT3
FT4 RAPID TEST KIT	25/50/100	VFIA-FT4
B-HCG TEST KIT	25/50/100	VFIA-BHCG



NUCLEIC ACID EXTRACTION KITS

Vision Biyoteknoloji offers both manual and automated DNA & RNA extraction solutions for different samples.

Our expert team utilizes the latest technology to generate simple, high-yield, time-saving and innovative nucleic acid purifying systems.

NUCLEIC ACID EXTRACTION KITS

PRODUCT	QUANTITY	REF
PURIFY AUTOMATED DNA EXTRACTION KIT	25/50/100	VIS100201
PURIFY AUTOMATED RNA EXTRACTION KIT	25/50/100	VIS100202
PURIFY DRIED BLOOD DNA EXTRACTION KIT S	25/50/100	VIS100203
PURIFY DRIED BLOOD RNA EXTRACTION KIT S	25/50/100	VIS100204
PURIFY BACTERIAL DNA EXTRACTION KIT M	25/50/100	VIS100231
PURIFY BACTERIAL RNA EXTRACTION KIT M	25/50/100	VIS100232
PURIFY BACTERIAL DNA EXTRACTION KIT S	25/50/100	VIS100233
PURIFY BACTERIAL RNA EXTRACTION KIT S	25/50/100	VIS100234
PURIFY GAITA DNA EXTRACTION KIT M	25/50/100	VIS100281
PURIFY GAITA RNA EXTRACTION KIT M	25/50/100	VIS100282
PURIFY GAITA DNA EXTRACTION KIT S	25/50/100	VIS100283
PURIFY GAITA RNA EXTRACTION KIT S	25/50/100	VIS100284
PURIFY PLANT DNA EXTRACTION KIT M	25/50/100	VIS100251
PURIFY PLANT RNA EXTRACTION KIT M	25/50/100	VIS100252
PURIFY PLANT DNA EXTRACTION KIT S	25/50/100	VIS100250
PURIFY PLANT RNA EXTRACTION KIT S	25/50/100	VIS100254
PURIFY MOLD-YEAST DNA EXTRACTION KIT M	25/50/100	VIS100241
PURIFY MOLD-YEAST RNA EXTRACTION KIT M	25/50/100	VIS100242
PURIFY MOLD-YEAST DNA EXTRACTION KIT S	25/50/100	VIS100243
PURIFY MOLD-YEAST RNA EXTRACTION KIT S	25/50/100	VIS100244
PURIFY OSTEOID DNA EXTRACTION KIT M	25/50/100	VIS100271
PURIFY OSTEOID RNA EXTRACTION KIT M	25/50/100	VIS100272
PURIFY OSTEOID DNA EXTRACTION KIT S	25/50/100	VIS100273
PURIFY OSTEOID RNA EXTRACTION KIT S	25/50/100	VIS100274
PURIFY TISSUE DNA EXTRACTION KIT M	25/50/100	VIS100261
PURIFY TISSUE RNA EXTRACTION KIT M	25/50/100	VIS100262
PURIFY TISSUE DNA EXTRACTION KIT S	25/50/100	VIS100263
PURIFY TISSUE RNA EXTRACTION KIT S	25/50/100	VIS100264
PURIFY VIRAL DNA EXTRACTION KIT M	25/50/100	VIS100221
PURIFY VIRAL RNA EXTRACTION KIT M	25/50/100	VIS100222
PURIFY VIRAL DNA EXTRACTION KIT S	25/50/100	VIS100223
PURIFY VIRAL RNA EXTRACTION KIT S	25/50/100	VIS100224
PURIFY WHOLE BLOOD DNA EXTRACTION KIT M	25/50/100	VIS100211
PURIFY WHOLE BLOOD RNA EXTRACTION KIT M	25/50/100	VIS100212
PURIFY WHOLE BLOOD DNA EXTRACTION KIT S	25/50/100	VIS100213
PURIFY WHOLE BLOOD RNA EXTRACTION KIT S	25/50/100	VIS100214
PURIFY SALIVA DNA EXTRACTION KIT S	25/50/100	VIS100215
PURIFY SALIVA RNA EXTRACTION KIT S	25/50/100	VIS100216

NUCLEIC ACID EXTRACTION KITS

PRODUCT	QUANTITY	REF
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T1	64	VIS100288
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T2	64	VIS100289
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T3	64	VIS100290
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T4	64	VIS100291
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T5	64	VIS100292
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T6	64	VIS100293
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T7	64	VIS100294
PURIFY AUTOMATED DNA&RNA EXTRACTION KIT T8	64	VIS100295
MANUEL VIRAL/ BACTERIAL EXTRACTION KIT FOR STOOL SAMPLES	25/50/100	VIS100296
MANUEL VIRAL/ BACTERIAL EXTRACTION KIT FOR URINE SAMPLES	25/50/100	VIS100297
MANUEL VIRAL/ BACTERIAL EXTRACTION KIT FOR BLOOD SAMPLES	25/50/100	VIS100298
MANUEL VIRAL/ BACTERIAL EXTRACTION KIT FOR SERUM SAMPLES	25/50/100	VIS100299



MICROBIOLOGY PANEL

Vision Biotechnology offers a wide range of ready-to-use sterile nutrient mediums for routine laboratory analysis and research purposes.

MICROBIOLOGY PANEL

PRODUCT	QUANTITY	REF
%5 SHEEP BLOOD AGAR	10 petri	VISB001
EMB AGAR	10 petri	VISB002
%5 SHEEP BLOOD AGAR/EMB AGAR	10 petri	VISB003
NUTRIENT AGAR	10 petri	VISB004
MUELLER HINTON AGAR	10 petri	VISB005
SABOURAUD DEXTROSE AGAR (SDA)	10 petri	VISB006
SALMONELLA SHIGELLA AGAR	10 petri	VISB007
SIM MEDIUM	10 tube	VISB008
SIMMONS CITRATE AGAR	10 tube	VISB009
TRIPLE SUGAR AGAR	10 petri	VISB010
UREA AGAR	10 tube	VISB011
CHOCOLATE AGAR	10 petri	VISB012
MAC CONKEY AGAR	10 petri	VISB013
MANNITOL SALT AGAR	10 petri	VISB014
LYSINE IRON AGAR	10 petri	VISB015
INDOLE BROTH	10 tube	VISB016
BILE ESCULIN AGAR	10 petri	VISB017
CROMAGAR VRE AGAR	10 petri	VISB018
CANDIDA CHROMOGENIC AGAR	10 petri	VISB019
%5 SHEEP BLOOD AGAR /EMB AGAR/SDA AGAR	10 petri	VISB020
%5 SHEEP BLOOD AGAR /EMB AGAR/CHOCOLATE AGAR	10 petri	VISB021
LÖWENSTEIN-JENSEN AGAR	10 petri	VISB022
MIO MEDIUM AGAR	10 petri	VISB023
SCHAEDLER AGAR	10 petri	VISB024
SDA/SDA/SDA AGAR	10 petri	VISB025
SS /SS /SS AGAR	10 petri	VISB026
MUELLER HILTON AGAR (120 MM)	10 petri	VISB031
MUELLER HILTON AGAR (150 MM)	10 petri	VISB032
CETRIMIDE AGAR	10 petri	VISB033
CLED AGAR	10 petri	VISB034
CHROMOGENIC COLIFORM AGAR	10 petri	VISB035
TRYPTIC SOY AGAR	10 petri	VISB036
SLANETZ & BARTLEY AGAR	10 petri	VISB037
BRILLANT GREEN AGAR	10 petri	VISB038

PRODUCT
QUANTITY
REF

COLUMBIA AGAR W/%5 SHEEP'S BLOOD	10 petri	VISB039
ENDO AGAR	10 petri	VISB041
EUGON LT 100 BROTH	1 bottle	VISB042
SORBITOL MAC CONKEY AGAR	10 petri	VISB043
MALT EXTRACT AGAR	10 petri	VISB044
MYCROPLASMA AGAR	10 petri	VISB045
NUTRIENT AGAR W/TTC	10 petri	VISB046
PLATE COUNT AGAR	10 petri	VISB047
POTATO DEXTROSE AGAR	10 petri	VISB048
RAMBACH AGAR	10 petri	VISB049
SDA W/ GENTAMICIN + CHLORAMPHENICOL	10 petri	VISB050
THAYER MARTIN AGAR	10 petri	VISB051
TSA W/NEUTRALISANT AGAR	10 petri	VISB052
VIOLET RED BILE AGAR	10 petri	VISB053
XLT4 AGAR	10 petri	VISB054
TSA W/PENASE AGAR	10 petri	VISB056
R2A AGAR	10 petri	VISB059
TSA W/BETALACTAMASE+ NEUTRALISANT (TRIPLE BAG)	10 petri	VISB060
OXFORD AGAR	10 petri	VISB061
ENTEROCOCCOSEL AGAR	10 petri	VISB062
D-COCCOSEL AGAR	10 petri	VISB063
TSC AGAR	10 petri	VISB064
DG18 AGAR	10 petri	VISB065
OGYE AGAR	10 petri	VISB066
TSA W/PENASE (TRIPLE BAG)	10 petri	VISB067
BIGGY AGAR	10 petri	VISB068
DNASE AGAR	10 petri	VISB069
CCDA AGAR	10 petri	VISB070
SDA W/CHLORAMPHENICOL	10 petri	VISB071
SCHAEDLER NEOMYCIN VANCOMYCIN AGAR W/%5 SHEEP BLOOD	10 petri	VISB072
MCTA (DOUBLE BAG)	10 petri	VISB073
MCTA W/BETALACTAMASE (DOUBLE BAG)	10 petri	VISB074
TSA W/PENASE+2 NEUTRALISANT (TRIPLE BAG)	10 petri	VISB075
TSA W/PENASE (DOUBLE BAG)	10 petri	VISB076

PRODUCT
QUANTITY
REF

SDA W/CHLORAMPHENOCOL + CYCLOHEXIMIDE)	10 petri	VISB077
D/E NEUTRALIZING AGAR (TRIPLE BAG)	10 petri	VISB078
%5 SHEEP BLOOD AGAR W/ BACITRACIN	10 petri	VISB079
COLUMBIA AGAR W/%5 SHEEP BLOOD / CHOCOLATE AGAR	10 petri	VISB080
CHROMAGAR VRE / CHROMAGAR MRSA	10 petri	VISB081
MANNITOL SALT AGAR / CETRIMIDE AGAR	10 petri	VISB082
XLD AGAR / M-ENDO LES AGAR (MELA)	10 petri	VISB083
SS/MACCONKEY AGAR	10 petri	VISB084
SDA W/ GENTAMYCIN + CHLORAMPHENIKOL / CHROMAGAR CANDIDA	10 petri	VISB085
%5 SHEEP BLOOD AGAR / MAC CONKEY AGAR w/o SALT + CRYSTAL VIOLET	10 petri	VISB086
CHROMOGENIC ORIENTATION AGAR/COLUMBIA CNA AGAR	10 petri	VISB087
MAC CONKEY AGAR/ COLUMBIA CNA AGAR	10 petri	VISB088
%5 SHEEP BLOOD AGAR / MUELLER HINTON AGAR	10 petri	VISB089
CHOCOLATE AGAR w/ISOVITALE X / SDA	10 petri	VISB090
MUELLER HINTON AGAR HORSE BLOOD	10 petri	VISB091
%5 SHEEP BLOOD AGAR w/%1,5 NaCl	10 petri	VISB092
LACTOSE TTC AGAR w/TERGITOL	10 petri	VISB093
SLANETZ BARTLEY (AZIDE) AGAR	10 petri	VISB094
BCYE AGAR	10 petri	VISB095
CETRIMIDE AGAR (60 MM PETRI)	10 petri	VISB096
VIOLET RED BILE GLUCOSE AGAR	10 petri	VISB098
BRAIN HEART INFUSION AGAR	10 petri	VISB099
SCHAEDLER AGAR w/%5 SHEEP BLOOD	10 petri	VISB100
ROGOSA AGAR	10 petri	VISB101
TSA NEUTRALISANT (TRIPLE BAG)	10 petri	VISB102
TSA DOUBLE BAG	10 petri	VISB103
SDA (DOUBLE BAG)	10 petri	VISB104
SDA CYCLOHEXIMIDE	10 petri	VISB105
TSA/PENASE + NEUTRALISANT	10 petri	VISB106
BAIRD PARKER AGAR	10 petri	VISB108
DRBC AGAR	10 petri	VISB110
EUGON AGAR	10 petri	VISB111
COLUMBIA AGAR	10 petri	VISB112
D/E NEUTRALIZING AGAR	10 petri	VISB113

PRODUCT	QUANTITY	REF
REINFORCED CLOSTRIDIAL AGAR	10 petri	VISB114
TBX AGAR	10 petri	VISB115
TRYPTIC SOY AGAR %5 HORSE BLOOD	10 petri	VISB116
VANCOMYCIN SCREEN AGAR	10 petri	VISB117
DESOXYCHOLATE CITRATE AGAR	10 petri	VISB118
LACTOBACILLUS MRS AGAR	10 petri	VISB119
%5 HUMAN BLOOD AGAR	10 petri	VISB120
ANAEROBIC AGAR	10 petri	VISB121
CHOCOLATE AGAR W/VBC	10 petri	VISB122
YERSINIA SELECTIVE (CIN) AGAR	10 petri	VISB124
THAYER-MARTIN AGAR (MODIFIED)	10 petri	VISB125
MYCOSEL AGAR	10 petri	VISB126
LAKED BLOOD AGAR	10 petri	VISB127
PHENYLETHYL ALCOHOL AGAR	10 petri	VISB128
SS AGAR / EMB AGAR	10 petri	VISB129
CHR.SALMONELLA AGAR/XLT4 AGAR	10 petri	VISB130
THIOGLYCOLLATE BROTH/7 ML	10 tube	VISB131
SELENITE CYSTINE BROTH/ 7 ML	10 tube	VISB132
TSB W/%0,05 SODIUM THIOGLYCOLLATE	10 petri	VISB133
ANTIBIOTIC AGAR NO 1	1 package	VISB134
ORNITHINE DECARBOXYLASE BROTH	10 tube	VISB135
UREA BROTH (CRISTENSEN)	10 tube	VISB136
COOKED MEAT MEDIUM	10 tube	VISB137
BUFFERED PEPTONE WATER (90 ML)	1 bottle	VISB138
TRYPTIC SOY BROTH	10 tube	VISB139
THIOGLYCOLLATE BROTH	10 tube	VISB140
MAC CONKEY BROTH	10 tube	VISB141
HALF FRASER BROTH	10 tube	VISB142
REINFORCED CLOSTRIDIAL MEDIUM	10 tube	VISB143
RAPPAPORT VASSILIADIS SALMONELLA ENRICHMENT BROTH	10 tube	VISB144
EE BROTH, MOSSEL	1 bottle	VISB145
BUFFERED SODIUM CHLORIDE SOLUTION	1 bottle	VISB146
FLUID THIOGLYCOLLATE MEDIUM	10 tube	VISB147
TSB STERILITY TEST BROTH	1 bottle	VISB148
FLUID A STERILITY TEST MEDIUM	1 bottle	VISB149

PRODUCT
QUANTITY
REF

ANTIBIOTIC BROTH NO 3	1 package	VISB150
TWEEN 80 0,1 %	1 bottle	VISB151
MYCOPLASMA BROTH BASE	1 package	VISB152
SEMISOLID NUTRIENT AGAR	10 petri	VISB153
MOTILITY INDOLE UREA MEDIUM	10 tube	VISB154
FLUID THIOGLYCOLLATE MEDIUM W/ BETA-LACTAMASE+TWEEN80	10 petri	VISB155
TRYPTIC SOY BROTH W/BETALACTAMASE	10 petri	VISB156
BUFFERED SODIUM CHLORIDE SOLUTION PH 7,0 W/ NEUTRALISANT	10 petri	VISB157
BUFFERED SODIUM CHLORIDE SOLUTION PH 7,0 W/ %0,3 LECITHIN + 2% TWEEN 80	1 bottle	VISB158
TRYPTIC SOY BROTH W/ %2 TWEEN 80 + %0,1 THYMIDINE	10 petri	VISB159
TSB + %4 TWEEN 20 + %0,5 LECITHIN	10 petri	VISB161
RVS BROTH	10 tube	VISB162
FLUID THIOGLYCOLLATE MEDIUM W/ TWEEN 80 (100 ML)	1 bottle	VISB163
FTM STERILITY TEST MEDIUM	10 tube	VISB165
TSB W-NEUTRALISANT	10 petri	VISB166
SDA W/CHLORAMPHENICOL	10 petri	VISB167
LOEFFLER MEDIUM	10 tube	VISB169
BISMUTH SULPHITE AGAR (WILSON- BLAIR)	10 petri	VISB170
HUGH LEIFSON MEDIUM	10 tube	VISB171
TODD HEWITT MEDIUM	10 tube	VISB172
MEAT EXTRACT W/ PEPTONE	10 petri	VISB173
PHENOL RED MEDIUM W/%1 GLUCOSE + %0,1 AGAR	10 petri	VISB174
BRAIN HEART INFUSION BROTH	10 tube	VISB175
BRAIN & BURKEY MEDIUM	10 tube	VISB176
CLARK-LUBS (MR-VP) BROTH	10 tube	VISB177
GIOLITTI-CANTONI BROTH	10 tube	VISB178
KLIGLER IRON AGAR	10 tube	VISB180
MUELLER HINTON BROTH (5 ML)	10 tube	VISB181
POTATO DEXTROSE AGAR (5 ML)	10 tube	VISB182
SABOURAUD DEXTROSE AGAR (5 ML)	10 tube	VISB183
SDA W/CHLORAMPHENICOL (5 ML)	10 tube	VISB184
SELENIT F BROTH (5 ML)	10 tube	VISB185
THIOGLYCOLLATE BROTH (5 ML)	10 tube	VISB186


PRODUCT
QUANTITY
REF

TRYPTIC SOY BROTH W/ %20 GLYCEROL (5 ML)	10 tube	VISB187
MAC CONKEY BROTH (5 ML)	10 tube	VISB188
BRILLANT GREEN BROTH (DURHAM TUBE) (5 ML)	10 tube	VISB189
RAPPAPORT VASSILIADIS (RVS) BROTH (9 ML)	10 tube	VISB190
NUTRIENT BROTH (5 ML)	10 tube	VISB191
GLUCOSE BROTH (5 ML)	10 tube	VISB192
SABOURAUD DEXTROSE BROTH (5 ML)	10 tube	VISB193
GN BROTH (5 ML)	10 tube	VISB194
COLUMBIA BLOOD BROTH (5 ML)	10 tube	VISB195
RVS BROTH (10 ML)	10 tube	VISB196
6,5 % NACL SOLUTION (5 ML)	10 tube	VISB197
LSTB W/MUG -DURHAM TUBE (9 ML)	10 tube	VISB198
MUELLER HINTON BROTH W/VANCOMYCIN + CEFTAZIDIME (5 ML)	10 tube	VISB199
MUELLER KAUFFMANN TETRATHIONATE NOVOBIOCINE BROTH (10 ML)	10 tube	VISB200
SKIM MILK MEDIUM (5 ML)	10 tube	VISB201
TSB W/NEUTRALISANT (GSTHLS) (10 ML)	10 tube	VISB202
THIOGLYCOLLATE BROTH (10 ML)	10 tube	VISB203
BILE ESCULIN BROTH (5 ML)	10 tube	VISB204
ENTEROCOCCOSEL BROTH W/ VANCOMYCIN (5 ML)	10 tube	VISB205
MOELLER DECARBOXYLASE BROTH W/ORNITHINE (5 ML)	10 tube	VISB206
FLUID THIOGLYCOLLATE MEDIUM (9 ML)	10 tube	VISB207
MOELLER DECARBOXYLASE BROTH (5 ML)	10 tube	VISB208
SDA W/GENTAMYCIN+ CYCLOHEXIMIDE (5 ML)	10 tube	VISB209
MUELLER HINTON BROTH (9 ML)	10 tube	VISB210
TRYPTIC SOY BROTH W/ BETA-LACTAMASE + TWEEN 80 (9 ML)	10 tube	VISB211
FLUID THIOGLYCOLLATE MEDIUM W/ BETA-LACTAMASE + TWEEN 80 (9 ML)	10 tube	VISB212
TRYPTIC SOY BROTH W/ TWEEN 80 (9 ML)	10 tube	VISB213
FLUID THIOGLYCOLLATE MEDIUM W/ TWEEN 80 (9 ML)	10 tube	VISB214
TRYPTIC SOY BROTH (10 ML)	10 tube	VISB215

PRODUCT**QUANTITY****REF**

SABOURAUD DEXTROSE AGAR (10 ML)	10 tube	VISB216
PSEUDOMONAS SELECTIVE AGAR	10 petri	VISB218
LETHEEN BROTH	10 tube	VISB219
TRYPTIC SOY BROTH	10 tube	VISB220
M-CP AGAR	10 petri	VISB222
PLATE COUNT AGAR	10 petri	VISB223
TBX CHROMOGENIC AGAR	10 petri	VISB224
VIOLET RED BILE (VRB) AGAR	10 petri	VISB225
DICHLORAN ROSE BENGAL CHLORAMPHENICOL	10 petri	VISB226
SALMONELLA CHROMOGENIC MEDIA	10 petri	VISB227
KF STERPTOCOCCUS AGAR	10 petri	VISB228
MANNITOL-EGG-YOLK-POLYMYXIN (MYP) AGAR	10 petri	VISB230
XYLOSE LYSINE DEOXYCHOLATE (XLD) AGAR	10 petri	VISB231
CHROMOGENIC LISTERIA AGAR	10 petri	VISB233
POTATO DEXTROSE AGAR (PDA)	10 petri	VISB235
NUTRIENT BROTH	10 tube	VISB236
BRILLIANT GREEN BILE (BGB) BROTH	10 tube	VISB237
COLIFORM AGAR	10 petri	VISB240
CHROMOGENIC CANDIDA AGAR	10 petri	VISB241
SHEEP BLOOD AGAR/UTI AGAR	10 petri	VISB242
CAMPLYBACTER AGAR	10 petri	VISB243
HEKTOEN ENTERIC AGAR	10 petri	VISB244
CHROMOGENIC KPC AGAR	10 petri	VISB245
CHROMOGENIC UTI AGAR	10 petri	VISB246
BLOOD CULTURE, BIPHASIC, PEDIATRIC	30 ml	VISB247
BLOOD CULTURE, BIPHASIC, ADULT	30 ml	VISB248

PRODUCT	QUANTITY	REF
DISTILLED WATER DD (PURE WATER)	500 ml / 1 lt / 2,5 lt	VISB249
IVD ISOTONIC SODIUM CHLORIDE SOLUTION	500 ml / 1 lt / 2,5 lt	VISB250
CH DILUENT	1,5 lt	VISB251
CH WASH	1,5 lt	VISB252



BLOOD GROUPING

Vision Blood grouping tests are gel column systems that provide easy and reliable blood grouping tests such as ABO grouping, Rh (D) identification, Cross match, An body screening, An body identification, Rh Subgroup determination (C, c, E, e and K).

GEL CARDS

PRODUCT	QUANTITY	REF
FORWARD ABO PATIENT	48	VIS-FP
FORWARD AB PLUS	48	VIS-FABP
FORWARD / REVERSE	48	VIS-FR
FORWARD / REVERSE DONORS	48	VIS-FRD
COMP / CROSS	48	VIS-CC
NEWBORN	48	VIS-NB
AHG	48	VIS-AHG
NEUTRAL / AHG	48	VIS-NAHG
RH PHENOTYPE	48	VIS-RP
RH PHENOTYPE PLUS	48	VIS-RPP
DAT DUO	48	VIS-DATDUO
DAT SCREEN	48	VIS-DAT
A1 GEL	48	VIS-A1
HGEL	48	VIS-H
NEUTRAL GEL	48	VIS-NEU
FORWARD / CDE	48	VIS-FCDE
IgG GEL	48	VIS-IGG

SOLUTIONS

PRODUCT	QUANTITY	REF
CELL REVERSE A1/B	2X10ML	VIS-RBC2
CELL SCREEN 1-2-3	3X5ML	VIS-RBCAS3
CELL SCREEN 1-2	2X5ML	VIS-RBCAS2
CELL SCREEN 1-2-3 PAP	3X5ML	VIS-RBCAS3P
CELL SCREEN 1-2 PAP	2X5ML	VIS-RBCAS2P
CELL COOMBS CONTROL	5 ML	VIS-CCC
BROMELIN	10 ML	VIS-BS10
LISS	500 ML	VIS-LISS500
QUALITY CONTROL CELL	4X5ML	VIS-CRBC
IgG SENSITIZED CELL	2 ML	VIS-DCRBC

GEL CARDS

PRODUCT

REF

GEL CARD CENTRIFUGE

VIS-JKS24

- 24 gel cards capacity
- Audible alarm
- Automatic lid opening system



GEL CARD READER

VIS-JKR1

- Computer connection
- Test result image
- Test results record
- Test results can be sent to the automation system or printer



GEL CARD INCUBATOR

VIS-JKI24

- 24 gel cards capacity
- Audible alarm
- Two independent heating sections





CONSUMABLES

CONSUMABLES

PRODUCT	QUANTITY	REF
PETRI DISH 90	700	VISPD90
PETRI DISH 90 2 SECTION	700	VISPD902
PETRI DISH 90 3 SECTION	700	VISPD903
PETRI DISH 90-ASEPTIC	700	VISPD90A
PETRI DISH 90 2 SECTION - ASEPTIC	700	VISPD902A
PETRI DISH 90 3 SECTION - ASEPTIC	414	VISPD903A
PETRI DISH 90-GAMMA STERILE	414	VISPD90G
PETRI DISH 90 2 SECTION - GAMMA STERILE	414	VISPD902G
PETRI DISH 90 3 SECTION- GAMMA STERILE	414	VISPD903G
100 ML URINE CONTAINER NON STERILE	1000	VISUC100
100 ML URINE CONTAINER ASEPTIC STERILE	400	VISUC100A
100 ML URINE CONTAINER GAMMA STERILE	200	VISUC100G
100 ML VACUUM URINE CONTAINER NON STERILE	500	VISVUC100
100 ML VACUUM URINE CONTAINER ASEPTIC STERILE	400	VISVUC100A
TEST TUBE(BLOOD) PP	5000	VISTTP1275
TEST TUBE (BLOOD) PS	5000	VISTTS1275
CONICAL TEST TUBE PP	2500	VISTTP16100
CONICAL TEST TUBE PS	2500	VISTTS16100
EPPENDORF TUBE	10000	VISET15
VACUUM URINE TEST TUBE	1000	VISVUT
STOOL CONTAINER WITH SPOON	2400	VISSCS
STOOL CONTAINER WITHOUT SPOON	1000	VISSCSW
TEST TUBE COVER	14000	VISTTC012
TEST TUBE COVER	10000	VISTTC016
TEST CASSETTE (1 HOLE)	3000	VISTT1
TEST CASSETTE (6 HOLE)	3000	VISTT6
TEST CASSETTE (12 HOLE)	3000	VISTT12
SLIDE BOX	1000	VISSB
SMEAR CYTOBRUSH	1000	VISSCB
SMEAR SAMPLING BRUSH	1000	VISSSB

VISION Q19-23 Real-Time PCR Fluorescent Quantitative System



VISION Q19-23 Real-Time PCR fluorescent quantitative system integrates a variety of advanced technologies, including using the thermoelectric

Peltier cooler module from United States, the maximum ramp rate is 6.5 °C/s, the available cycles is more than 1 million times, and it is compatible with win10 operating system. It is specially designed to meet the needs of rapid nucleic acid detection in clinical, food, environmental and research applications.

Options of 1-4 fluorescence channels and 8-16 samples can be flexibly customized, and can be freely matched with experiments to achieve realtime and high-sensitivity fluorescence measurement.

Flexible Selection of Fluorescent Channels

Single Channel (FAM)

Dual Channel (FAM/VIC)

Four Channels (FAM/VIC/ROX/CY5)

MODEL	VISION Q19-23 8	VISION Q19-23 16	
Throughput	8 samples	16 samples	
Consumable	0.2ml thin-walled clear tube, 8x0.2ml thin-walled clear tube		
Reaction Volume	20-100 μ l		
Temperature Control Technology	Thermoelectric Peltier cooler module		
Heating Temperature Range	0-100°C(Accuracy: 0.1°C)		
Ramp Rate	Maximum ramp rate: 6.5°C/s Average ramp rate: not lower than 4.5°C/s		
Temperature Control Accuracy	0.1°C		
Temperature Uniformity	\pm 0.2°C		
Temperature Accuracy	\pm 0.2°C		
Temperature Range Of Hot Lid	30-110°C (Adjustable, default 105°C)		
Temperature Control Mode	Block/Tube		
Fluorescence Intensity Repeatability	CV \leq 1%		
Fluorescence Intensity Precision	CV \leq 3%		
Fluorescence Channel	4 channels (FAM/VIC/ROX/CY5)	2 channels(FAM/VIC)	1 channel(FAM)
Fluorescence Excitation Light Wavelength	F1: 465nm, F2: 525nm F3: 571nm, F4: 628nm	F1: 465nm F2: 525nm	F1: 465nm
Fluorescence Detection Wavelength	F1: 510nm, F2: 564nm F3: 612nm, F4: 692nm	F1: 510nm F2: 564nm、	F1: 510nm
Excitation Light Source	High brightness LED, long service life, no maintenance		
Detector	Photodiode (PD)		
Dynamic Range	1-10 ¹⁰ Copies		
Detection Sensitivity	Can detect single copy genes		
Sample Detection Repeatability	CV \leq 3%		
Linearity	r \geq 0.990		
Software Analysis Function	Quantitative/Qualitative, Melting Curve, Genotyping		
Print	Supports direct printing of reports		
Control Method	External computer control, Win10 operating system		
Interface	USB2.0		
Dimension	width x depth x height=220×200×150mm		
Weight	3.5kg		
Input Voltage	DC15V		
Output Power	240W		

VIS-MAG EXTRACTION SYSTEM 32



Specification

- Method: Magnetic bead
- Sample type: Whole blood, serum, plasma, throat swab, puncture fluid, etc.
- Throughput: 32 samples
- Extraction volume: 20-1000 μ L
- Magnetic beads recovery rate: 95%
- Extraction time: 9-18minutes

Efficient

Unique extraction mode , easy and high efficiency of batch operation with short extraction mode (9-18 minutes).

Precision

Precise constant temperature control system and precise design of magnetic beads to assure the high purity extraction.

Safety

Extraction process is completed in a closed environment and disinfected by ultraviolet lamp to avoid pollution.

Convenient

Friendly-use so ware and different extraction methods can be set according to the sample type.



VISION BIYTEKNOLOJI LTD.



visionbiotechnology



/vision-biotechnology/



+90 2167840298



info@visionbiotechnology.com



<https://www.visionbiotechnology.com/>



Yukarı Dudullu Mah. Alemdağ Cad. Hıra İş Hanı
Blok No:742a Ümraniye/ İstanbul

©2025 Vision Biyoteknoloji. All rights reserved. All trademarks are the property of Vision Biyoteknoloji and its subsidiaries. Any photos displayed are for illustrative purposes only. Specification and terms are subject to change. VSN.2024.02



INFINITY IVD